

*No. 1 choice  
of hospitals  
& mothers*

# Optimise Hospital Processes with Ready-to-Use Products



- I **Ready-to-Use:** No cleaning needed prior to first use
- I **Safe:** Always new and therefore hygienic
- I **Simple:** Disposable after use
- I **Convenient:** More time available to spend with babies in your care



# Improve your Hospital Process – with our simple and Ready-to-Use Products.

Our breastfeeding products are the result of more than 30 years of evidence-based research. Thanks to this accumulated knowledge, Medela is the partner for complete solutions in the hospital sector. Our pump sets, breastmilk bottles and teats help to get precious breastmilk safely to where it is needed: in the baby's tummy. After all, we all want the same thing: happy and healthy mothers and babies.

With its Ready-to-Use product range, Medela greatly simplifies the daily routine involved in breastfeeding. All these products can be used straight away without cleaning prior to first use. This supports hygiene and is safe to use. Cost savings are likely due to the elimination of the time, cost and energy-intensive sterilisation/ disinfection process.



**Breastfeeding Solution Circle**

## Safe and gentle: Medela pump sets.


Our pump sets are designed for optimised and convenient handling. Every mother gets a new personal set to minimise the risk of breastmilk contamination. Since the One-Day Pump Sets can be used for up to eight pumping sessions within 24 hours, one set a day is enough. During this time the mother takes responsibility for her pump set (e.g. cleaning, storing). The simplicity of the design makes it easy for her to clean and handle it. This saves you costs and time which you can devote to your other duties and attending to other patients. To satisfy the mothers' individual needs we offer different breastshield sizes to ensure maximum comfort and efficiency while pumping. Moreover, the appropriate bottle size can be chosen according to the personal milk output.

## Handy and hygienic: Medela breastmilk bottles.

Medela Disposable Bottles are used to collect, store, portion and feed breastmilk. In a milk bank, the bottles can also be used for pasteurisation. We provide different sizes of bottles: 35ml, 80ml and 150ml.

The 35ml Colostrum Container, for example, is designed with healthcare professionals and mothers in mind. The curved bottom minimises the loss of colostrum or breastmilk during transfer to a syringe. The smaller size helps to keep mothers motivated to pump by setting realistic expectations of initial milk production.


The breastmilk bottles are graduated and enable easily to check the quantity of breastmilk that has been expressed. In addition, you can note the baby's name, the quantity of breastmilk expressed, as well as the time and date of expression on the additional available label.



**One-Day Pump Sets**



**Disposable Breastmilk Bottles**



**Label for Disposable Breastmilk Bottles**

### Individually fit: Medela teats.

Our practical, individually packed disposable teats ensure that it is simple to feed a baby with breastmilk. Two different milk flow rates accommodate the different needs of premature and newborn babies. As an alternative to teat feeding, Medela offers a wide range of special feeding devices, for example devices to support babies with cleft lip/palate abnormalities.



Medela breastmilk bottles and all our products that come in contact with breastmilk are and always have been BPA-free.



## Medela's Reusable Products

For hospitals with traditional cleaning and sterilising practices.

### Pump Sets

- I Designed for use by multiple mothers
- I Autoclavable
- I Single and Double Pump Sets available
- I Various breastshield sizes available (S, M, L, XL, XXL)



Symphony Pump Set



Basic Kit

### Breastmilk Bottles

- I Designed for use by multiple mothers
- I Autoclavable
- I Sizes 80ml and 150ml available



80 ml



150 ml

### Teats

- I Autoclavable
- I Slow and Medium Flow Teat available
- I Available in bulk packs (25 pcs.)



Slow Flow



Medium Flow

Ready-to-Use Products			
<b>Pump Sets</b>			
Article-Number	Product	Materials <sup>1)</sup>	Qty
008.0176	Symphony One-Day Pump Set with PersonalFit Size M (24mm)	PP, TPE	54
008.0177	Symphony One-Day Pump Set with PersonalFit Size L (27mm)	PP, TPE	45
008.0178	Symphony One-Day Pump Set with PersonalFit Size XL (30mm)	PP, TPE	45
008.0179	Symphony and Lactina One-Day Pump Set with PersonalFit Size M (24mm)	PP, TPE	72
008.0180	Symphony and Lactina One-Day Pump Set with PersonalFit Size L (27mm)	PP, TPE	60
008.0181	Symphony and Lactina One-Day Pump Set with PersonalFit Size XL (30mm)	PP, TPE	60
<b>Breastmilk Bottles</b>			
200.2110	35ml Colostrum Container	PP, PBT	100
200.2957	80ml Disposable Breastmilk Bottles	PP, HDPE	480
200.2956	150ml Disposable Breastmilk Bottles	PP, HDPE	320
<b>Teats</b>			
200.3200	Disposable Hospital Teat, Standard, Medium Flow	PP, TPE	100
200.3201	Disposable Hospital Teat, Premature, Slow Flow	PP, TPE	100
<b>Reusable Products</b>			
<b>Pump Sets</b>			
800.0557	Symphony Hospital Multibox with PersonalFit Size M (24mm)	PP, Silicone	20
200.1141	Symphony Hospital Multibox with PersonalFit Size L (27mm)	PP, Silicone	20
200.1143	Symphony Hospital Multibox with PersonalFit Size XL (30mm)	PP, Silicone	20
800.0605	Basic Kit (Breastshield M (24mm), valve, bottle, lid)	PP, Silicone	25
<b>Pump Set Accessories</b>			
200.3887	Label for Disposable Breastmilk Bottles	-	5x2'500
008.0182	Label Dispenser	-	1
800.0608	Symphony Membrane Cap and Protective Membrane	-	25
800.0609	Lactina Piston/Cylinder incl. Cylinder Holder	-	5
800.0659	Silicone Tubing for Symphony and Lactina	-	25
200.0571	PersonalFit Breastshield Size S (21mm)	PP	25
800.0668	PersonalFit Breastshield Size M (24mm)	PP	25
800.0711	PersonalFit Breastshield Size L (27mm)	PP	25
800.0712	PersonalFit Breastshield Size XL (30mm)	PP	25
800.0854	PersonalFit Breastshield Size XXL (36mm)	PP	25
<b>Breastmilk Bottles</b>			
810.0137	80ml Reusable Breastmilk Bottle	PP	1
800.0625	150ml Reusable Breastmilk Bottle	PP	25
800.0627	Lid	PP	25
800.0628	Disc	PP	25
<b>Teats</b>			
008.0119	Reusable Hospital Teat, Slow Flow	Silicone	25
008.0120	Reusable Hospital Teat, Medium Flow	Silicone	25

Medela also offers Hospital Grade Breastpumps: Symphony and Lactina. Find out more at [www.medela.com](http://www.medela.com) or contact your local Medela representative.  
<sup>1)</sup> in contact with breastmilk